

尺寸：148x210mm
 材质：128g双铜
 工艺：双面印刷+四色

feellife
 The third way of dosing



Air Plus III
 Portable Variable Frequency Nebulizer
 One nebulization is worth three times

Free switch between micro-wave mode and uniform-rate mode
 It can be suitable for patients at different ages

Micro-wave mode
 The nebulization rate changes circularly according to human respiration, bringing a more comfortable nebulization experience

Variable frequency mode
 The nebulization rate can be adjusted manually to improve the absorption


Self-cleaning technology
 To prevent the micro mesh blockage

Less residual liquid and low noise
 Making it a pleasant time during nebulization



Technical parameters


Model: Air Plus III
 Power Supply: Internal power supply
 Power Consumption: <4.0w
 Nebulization Rate: 0.15-0.90ml/min
 (Multi-frequency adjustable)
 MMAD: 2.5±30%µm
 Medication Cup Capacity: 10ml
 Device Dimensions:
 Air Plus III(43(L) * 55(W) * 99(H)mm)



Deposition site: reach to bronchia and alveoli


Authoritative test report of proportion of particle size

The medium position particle size of the nebulizers in the temperature of 25±3°C and the humidity of 60±5%RH with the 0.9% physiological saline. Under these conditions of the particle size distribution curve s as follows:



Notes: The horizontal axis is the particle size which hews Logarithmic distribution. The vertical axis on the left is the cumulative percentage of the volume, corresponding to the upward trend of the graph. The vertical axis on the right is the volume percentage of a certain interval, corresponding to the histogram or undulating curve.

* Application: To make the liquid medication nebulized so that the patients can inhale directly.



FEELLIFE HEALTH INC.
 Factory: Building A, NO. 9 Furong Road, Songgang Street,
 Bao'an District, Shenzhen 518104
 Website: www.feellife.com